

## Chem 131A: Folin-Ciocalteu Analysis for Protein

This method is quite sensitive. Samples containing as little as 5  $\mu\text{g}$  of protein can be readily analyzed.

The principle depends on the reaction of the phenolic (R) group of tyrosine with the colorless phosphomolybdate ions in the presence of  $\text{Cu}^{2+}$  ions. The phenol is oxidized and the complex phosphomolybdate is reduced to a blue chromogen that can be estimated spectrophotometrically.

**Important!! Steps 1 and 2 must be done at the same time.**

1) Standard Tyrosine Curve: To obtain data for the curve, set up the following tubes:

Tube	1	2	3	4	5	6	UNK	UNK
Standard tyrosine (mL)	0	0.10	0.20	0.30	0.40	0.50	0.30	0.50
Water (mL)	1.20	1.10	1.00	0.90	0.80	0.70	0.90	0.70
Absorbance								

Standard tyrosine solution: 100  $\mu\text{g}$  tyrosine/mL

2) Tyrosine Unknown.

You will receive an unknown containing tyrosine plus an inert compound. The percentage of tyrosine varies from 10 to 90 %. Weigh out an exact amount (calculate) of unknown in a beaker on an analytical balance ( $\sim 1\text{mg}$ ). Dissolve, then dilute to volume in a 10 mL volumetric flask. Set up two tubes as in the table above. Accurately determine the % tyrosine (3 sig. places) in the unknown.

When all 8 tubes have been prepared (this can be shared between two people):

1) Prepare 100 mL of fresh alkaline copper reagent by mixing together in the following order: 1 mL of 1%  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ , 1 mL of 2% sodium tartrate and 98 mL of 2%  $\text{Na}_2\text{CO}_3$  in 0.1 M NaOH.

2) Add and mix immediately 6 mL of fresh alkaline copper solution to each tube. Allow all the tubes to stand for 10 minutes at room temperature. Add and mix immediately 0.3 mL of Folin-Ciocalteu reagent to each tube. Allow 30 minutes for full color development. Measure absorbance at 500 nm. (Adapted from Clark and Switzer, EXPERIMENTAL BIOCHEMISTRY, 2<sup>nd</sup> edition, p.12, 1977, Freeman.)

3) Plot a standard curve: absorbance vs  $\mu\text{g}$  of tyrosine. Use the standard curve to determine  $\mu\text{g}$  of tyrosine in your unknown.

**Important!! Steps 3 and 4 must be done at the same time.**

3) Standard protein curve. To obtain data for this curve, set up the following tubes:

Tube	1	2	3	4	5	6	UNK	UNK
Standard Egg Albumin (mL)	0	0.10	0.20	0.30	0.40	0.50	0.30	0.50
Water (mL)	1.20	1.10	1.00	0.90	0.80	0.70	0.90	0.70
Absorbance								

Standard Egg Albumin Solution: 400  $\mu\text{g}$  protein/mL

4) You will be given a solution containing egg albumin at a concentration ranging from 70 to 380  $\mu\text{g}/\text{mL}$ . Set up two tubes as indicated in Table 2. Determine and report the concentration of the unknown protein.

When all 8 tubes have been prepared (this can be shared by two people):

1) Add and mix in the same reagents to each of these tubes as described for tyrosine determination (Items 1 and 2 under A).

2) Plot a standard curve: absorbance vs  $\mu\text{g}$  of protein. Use the standard curve to determine  $\mu\text{g}$  of protein in your unknown.